smiths medical
portex®

Spinal Needle Systems

- en** Instructions for Use
- fr** Mode d'emploi
- de** Gebrauchsanleitung
- it** Istruzioni per l'uso
- es** Instrucciones de uso
- pt** Instruções de utilização
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- da** Brugsanvisning
- no** Bruksanvisning
- fi** Käyttöohjeet
- el** Οδηγίες χρήσης
- cs** Návod k použití
- hu** Használati utasítás
- pl** Instrukcja stosowania
- et** Kasutusjuhiseid
- lv** Lietošanas instrukciju
- ro** Instrucțiuni de utilizare

 **Manufacturer:**
Smiths Medical ASD, Inc.
10 Bowman Drive
Keene, NH 03431 USA
Tel: 1 800 258 5361 (US/CA)
Tel: 1 614 210 7300

EC REP European Representative:
Smiths Medical International, Ltd.
1500 Eureka Park, Lower Pemberton
Ashford, Kent, TN25 4BF, UK
Tel: +44 (0)1233 722100
www.smiths-medical.com

Rx CE
ONLY 0086

2018-03

1 0019000-001

Portex® Spinal Needle Systems

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:

100/496/022	<i>Spinal Needle Set, Lancet Point Spinal Needle 22G</i>
100/496/024	<i>Spinal Needle Set, Lancet Point Spinal Needle 24G</i>
100/496/025	<i>Spinal Needle Set, Lancet Point Spinal Needle 25G</i>
100/496/026	<i>Spinal Needle Set, Lancet Point Spinal Needle 26G</i>
100/496/027	<i>Spinal Needle Set, Lancet Point Spinal Needle 27G</i>
100/496/122	<i>Spinal Needle Set, Pencil Point Spinal Needle, Size 22G</i>
100/496/124	<i>Spinal Needle Set, Pencil Point Spinal Needle, Size 24G</i>
100/496/125	<i>Spinal Needle Set, Pencil Point Spinal Needle, Size 25G</i>
100/496/126	<i>Spinal Needle Set, Pencil Point Spinal Needle, Size 26G</i>
100/496/127	<i>Spinal Needle Set, Pencil Point Spinal Needle, Size 27G</i>
100/492/815	<i>Spinal Needle Set, Extra-length Pencil Point Spinal Needle, Size 25G</i>
100/492/816	<i>Spinal Needle Set, Extra-length Pencil Point Spinal Needle, Size 26G</i>
100/492/817	<i>Spinal Needle Set, Extra-length Pencil Point Spinal Needle, Size 27G</i>

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

DESCRIPTION

A range of sterile, single use packs to perform a single injection of local anaesthetic or other drugs into the subarachnoid space for pain relief.

Spinal Midi-Trays contain in addition to the above:

- 1 Straw Filter (5µm)
- 1 Skin Wheel Needle 25G x 5/8" (0.5mm x 16mm)
- 1 Infiltration Needle 22G x 1 1/2" (0.7mm x 38mm)
- 1 3ml Plastic Syringe
- 1 5ml Plastic Syringe

INSTRUCTIONS FOR USE

- 1 Check that the spinal needle will pass freely through the introducer needle.
- 2 Perform aseptic preparation of the skin.
- 3 Infiltrate local anaesthetic into the skin, subcutaneous and deeper tissues.
- 4 Insert the introducer needle, avoiding excessive insertion which might pierce the dura. Check that the stylet of the spinal needle is correctly fitted. Insert the spinal needle through the introducer needle, and advance it through the resistances of

first the ligamentum flavum, and then the dura mater. After piercing the dura, the needle should be rotated so that the lateral eye (pencil point) is less likely to be occluded by a nerve root.

- 5 Withdraw the stylet from the spinal needle and wait for the cerebrospinal fluid (CSF) to appear in the needle hub. Attach the syringe containing local anaesthetic and aspirate gently to confirm free flow of CSF. If CSF does not appear in the needle hub, replace the stylet and rotate the needle 90° or advance or withdraw the needle until CSF appears and confirm free flow of CSF in the syringe.
- 6 Inject the local anaesthetic, withdraw the spinal and introducer needles, and cover the puncture site with a suitable dressing.

PRECAUTIONS

- 1 Ensure that the bevel on lancet point needles is facing laterally during puncture of the dura.
- 2 Always advance or re-position the spinal needle with the stylet in place.
- 3 If there is any resistance to injecting the local anaesthetic through the spinal needle, or pain on injection, stop immediately and proceed in accordance with currently accepted medical techniques.

WARNINGS

- 1 **Never use excessive force to advance the spinal needle when bone is felt. Continued advancement may bend the tip of the spinal needle. If undue resistance is felt, carefully remove the spinal needle and discard.**
- 2 **If resistance is encountered when removing or re-positioning the spinal needle, the spinal needle and introducer needle must be withdrawn together and the procedure repeated using new needles.**
- 3 **Always use spinal needles with an introducer or Tuohy needle to reduce the risk of the spinal needle bending. The tip of a bent spinal needle may break inside the patient when attempt is made to withdraw it. A bent spinal needle, or the surgical intervention required to remove a broken needle tip, will cause unnecessary trauma to the patient.**

Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2018 Smiths Medical. All rights reserved.

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en**Single use.**

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

Contents sterile unless package has been opened or damaged.

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Portex, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Portex, znak towarowy Portex oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.