DASH 6 SYRINGES

GB DASH 6™ ISO 80369-6 NRFit® Compliant Syringe

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GB

INTENDED USE

The DASH 6™ NRFit® Lock Syringe and DASH 6™ NRFit® Slip Syringes are intended to be used with ISO 80369-6 NRFit® neuraxial compliant devices for the administration of neuraxial medication or anaesthetic.

INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial device and is used to administrate neuraxial medication and anaesthetic to a patient.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for intravenous or enteral use.

WARNINGS & PRECAUTIONS

Do not use if packaging is damaged.

Do not use if there are no dose markings on the syringe.

Do no re-use or reprocess the device as this may lead to leakage or infection.

Do not use if the tip is visibly blocked or the device is damaged.

The device is single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device.

Verify that the neuraxial device being used is compatible with DASH 6™ ISO 80369-6 NRFit® to reduce the risk of misconnections.

Fully depress the plunger before filling with neuraxial medication/anaesthetic. After filling, ensure that the syringe moat is free of excess neuraxial medication/anaesthetic and that there are no air bubbles in the syringe. If there is excess neuraxial medication/anaesthetic or there are any air bubbles, remove them by flicking the side of the syringe.

Attach the ISO 80369-6 NRFit® compliant syringe into the compatible medication port.

Once this is secure, gently push the syringe plunger to transfer the medication/anaesthetic into the neuraxial site.

After completely expelling the medication/anaesthetic from the syringe, detach the syringe and dispose of immediately as per hospital protocol.

It is recommended to use a compatible filter in conjunction with the neuraxial device, or a drawing up device with an integrated filter, to prevent entrainment of any foreign particulate matter when drawing up fluid into the device.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

GRADUATION ACCURACY

Graduation accuracy detailed below complies with the requirements of ISO 7886-1.

Syringe volume (Nominal Capacity)	Capacity tolerance (greater than or half-filled syringe)	Capacity tolerance (less than half-filled syringe)
1ml	± 5% of expelled volume	± 1.5% of nominal syringe capacity x ± 2% of expelled volume
3ml	± 5% of expelled volume	± 1.5% of nominal syringe capacity x ± 2% of expelled volume
5ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
10ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
20ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume
60ml	± 4% of expelled volume	± 1.5% of nominal syringe capacity x ± 1% of expelled volume

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