smiths medical

portex®

Epidural Tuohy Needles

- en** Instructions for use
- fr** Mode d'emploi
- de** Gebrauchsanweisung
- it** Istruzioni per l'uso
- es** Instrucciones de uso
- pt** Instruções de utilização
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- da** Brugsanvisning
- no** Bruksanvisning
- fi** Käyttöohjeet
- el** Οδηγίες χρήσης
- cs** Návod k použití
- hu** Használati utasítás
- pl** Instrukcja użycia

 Manufacturer:
Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442 USA
Tel: 1 800 258 5361 (US/CA)
Tel: +1 614 210 7300

 European Representative:
Smiths Medical Czech Republic a.s.
Olomoucká 306, Hranice 1 - Město,
753 01 Hranice, Czech Republic
Tel: +44 (0) 1233 722100
www.smiths-medical.com

Rx **CE**
ONLY **2797**

2019-02

10019003-003

Portex® Epidural Tuohy Needles

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® products:

- 100/395/160 *Tuohy Needle, Size 16G (O.D. 1.6mm, I.D. 1.15mm, length 80mm)*
- 100/395/560 *Tuohy Needle, Extra Long, Size 16G (O.D. 1.6mm, I.D. 1.15mm, length 110mm)*
- 100/395/180 *Tuohy Needle, Size 18G (O.D. 1.3mm, I.D. 1.0mm, length 80mm)*
- 100/395/580 *Tuohy Needle, Extra Long, Size 18G (O.D. 1.3mm, I.D. 1.0mm, length 110mm)*

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

DESCRIPTION

A range of sterile, single use Tuohy needles for insertion into the epidural space for either single shot doses of local anaesthetic or other drugs or as an introducer for epidural catheters. The Tuohy needles comply with BS 6196, are graduated at 10mm intervals to enable the depth of needle insertion to be determined. The Tuohy needles are supplied with snap lock wing which allows the needle to be used as a winged or non-winged needle. The needle hub is colour coded for ease of identification.

INSTRUCTIONS FOR USE

1. If the Tuohy needle is being used to introduce an epidural catheter, check the catheter will pass through the needle.
2. Attach snap-lock wings (if required) ensuring the orientation lugs on the wings and needle hub are aligned.
3. Perform aseptic preparation of the skin.
4. Infiltrate local anaesthetic into the skin, subcutaneous and deeper tissues.
5. Identify the epidural space in accordance with currently accepted medical techniques.
6. Inject drugs of choice, remove needle and cover puncture site with suitable dressing.

PRECAUTIONS

1. If blood or cerebrospinal fluid (CSF) is observed during the placement of the Tuohy needle in the epidural space, the needle has been placed incorrectly. Stop immediately and proceed in accordance with currently accepted medical techniques.
2. It is recommended that a filter be used when injecting drugs into the epidural space by either Tuohy needle or epidural catheter to prevent particles or organisms being injected into the epidural space.

WARNINGS

Never pull back the catheter through the epidural Tuohy needle as this can result in the catheter being cut and left in the epidural space. If catheter insertion difficulties are encountered, the epidural Tuohy needle and catheter should be removed carefully together as a single unit, and the procedure repeated.

Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical. All rights reserved.

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Portex, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Portex, znak towarowy Portex oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol ® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.